

INTEOS[®] Olecranon Plate 3.0 & 3.5


Innovative osteosynthesis system with Olecranon plates



Advantages of the INTEOS® 3.0 & 3.5 Olecranon Plates system

One of the strengths of Hofer Medical Solutions INTEOS® Olecranon 3.0 & 3.5 plates is the particularly simple implantation technique: the standards used in this case are angular and polyaxial stability techniques.

- + Correct anatomical drawing
- + Angular and polyaxial stability
- + The hemispherical screw heads are concealed in the plate
- + Automatic compression holes to ensure the necessary compression
- + Numerous holes for K. threads to guarantee temporary fixation of the fragments
- + Induction to the consolidation process
- + Various plate lengths (3-5-7-9-11 holes)
- + Two different models - with and without loops - to improve and achieve maximum fracture stability



INTEOS® Olecranon Plate 3.0 & 3.5 as av



Titanium	Number of Holes	Length	Side	PU
730-135-000-012	9 Head / 3 Shaft	89 mm	left	1
730-135-000-014	9 Head / 5 Shaft	114 mm	left	1
730-135-000-016	9 Head / 7 Shaft	130 mm	left	1
730-135-000-018	9 Head / 9 Shaft	146 mm	left	1
730-135-000-020	9 Head / 11 Shaft	161 mm	left	1
730-135-001-012	9 Head / 3 Shaft	89 mm	right	1
730-135-001-014	9 Head / 5 Shaft	114 mm	right	1
730-135-001-016	9 Head / 7 Shaft	130 mm	right	1
730-135-001-018	9 Head / 9 Shaft	146 mm	right	1
730-135-001-020	9 Head / 11 Shaft	161 mm	right	1

Additional options

S	xxx-xxx-xxx-xxx- S	Sterile	
T	xxx-xxx-xxx-xxx- T	Trial	Do not implant!

INTEOS® Olecranon Plate 3.0 & 3.5 with loop as av



Titanium	Number of Holes	Length	Side	PU
730-135-100-015	9 Head / 6 Shaft	89 mm	left	1
730-135-100-017	9 Head / 8 Shaft	114 mm	left	1
730-135-100-019	9 Head / 10 Shaft	130 mm	left	1
730-135-100-021	9 Head / 12 Shaft	146 mm	left	1
730-135-100-023	9 Head / 14 Shaft	161 mm	left	1
730-135-101-015	9 Head / 6 Shaft	89 mm	right	1
730-135-101-017	9 Head / 8 Shaft	114 mm	right	1
730-135-101-019	9 Head / 10 Shaft	130 mm	right	1
730-135-101-021	9 Head / 12 Shaft	146 mm	right	1
730-135-101-023	9 Head / 14 Shaft	161 mm	right	1

Additional options

S	xxx-xxx-xxx-xxx- S	Sterile	
T	xxx-xxx-xxx-xxx- T	Trial	Do not implant!





AUSTRIA

HOFER GMBH & CO KG
Jahnstrasse 10-12
A 8280 Fürstenfeld
TEL: +43(0)3382/53388
FAX: +43(0)3382/53093
office@hofer-medical.com
www.hofer-medical.com
FN21826y
UID ATU30764704



GERMANY

HOFER-MEDICAL GMBH
Ruhrtalstraße 52-60
D 45239 Essen
TEL: 0800-700 80 22 (innerhalb BRD)
FAX: +49 201 56599 1069
office@hofer-medical.de
www.hofer-medical.com
HRB24645
UID DE290596924



ITALY

HOFER MEDICAL ITALIA S.R.L.
Via della Rena, 26
I - 39100 Bolzano

Tel: +39 393 8950481
FAX: +39 06 81151148
office@hofer-medical.it
www.hofer-medical.com
RM-1454597
UID IT13534721009



POLAND

HOFER GMBH & CO KG SPÓŁKA
KOMANDYTOWA ODDZIAŁ W POLSCE
Radkowska 9
PL 57-402 Nowa Ruda

Tel.Kom.:+48603224112
Fax: +48 748712273
office@hofer-medical.pl
www.hofer-medical.com
REGON 141527293
UID PL1070011771

creating partnership

For us, a partnership is the result of a long and happy relationship in all areas of our work . Reliable bone healing for patients, optimal handling of the systems for surgeons and surgical personnel, as well as the simple preparation of the instruments, constantly strengthen this partnership.

For us, "creating" means more than just finding solutions in the form of innovative products.

It is the result of high standards, constant development, innovative products and excellent service for patients, surgeons and surgical personnel around the clock.

Please do not hesitate to contact us if you have questions about our company, our employees or our production methods.

Your Hofer Medical SolutionsTeam

vCard

